

Priming and Cleaning Your Maxair™ Autohaler™

Patient Education Guide

Priming or Preparing Instructions for Your Maxair™ Autohaler™

Priming your Autohaler™ should be done before using it for the first time. If it has been more than 48 hours since you last used your Autohaler™, you will need to prime it again. Ask your health-care provider or pharmacist how many puffs of medicine your Autohaler™ has when it is full. You need to keep track of how many puffs of medicine you take everyday, so you can have your Autohaler™ refilled before you run out of medicine.



1 Remove the mouthpiece cover.



2 Point mouthpiece away from you.
Push lever up.



3 Push white slide on bottom of the
mouthpiece. Arrows on the bottom
will show you the direction to push the
slide.



4 Push lever back down then repeat
steps 2-3.

The American College of Chest Physicians is the leading resource for the improvement of cardiopulmonary health and critical care worldwide. Its mission is to promote the prevention and treatment of diseases of the chest through leadership, education, research, and communication.

This publication's content contains general information, is not intended to be and is not complete, is not medical advice, and does not replace professional medical care and physician advice, which always should be sought for any specific condition. The American College of Chest Physicians and its officers, regents, executive committee, members, and employees specifically disclaim all responsibility for any liability, damages (actual or consequential), loss, or risk, personal or otherwise, based on any legal theory whatsoever, alleged to have been incurred as a result, directly or indirectly, of the use of any of the material herein.

Priming and Cleaning Your Maxair™ Autohaler™

Page 2

AMERICAN COLLEGE OF
CHEST
P H Y S I C I A N S ®

American College of Chest Physicians
3300 Dundee Road, Northbrook, IL 60062
(847) 498-1400 phone
(847) 498-5460 fax
www.chestnet.org



Cleaning Your Maxair™ Autohaler™



1 Remove the mouthpiece cover.



2 Wipe the mouthpiece with a clean dry cloth.



3 Turn the Autohaler™ upside down. Tap back of the Autohaler™. Flap should come down so that the spray hole can be seen.



4 Clean flap with a dry cotton swab.



5 Put mouthpiece cover back on. Make sure lever is down. Clean the Autohaler™ once a week or more often if Autohaler™ is visibly dirty.

The American College of Chest Physicians is the leading resource for the improvement of cardiopulmonary health and critical care worldwide. Its mission is to promote the prevention and treatment of diseases of the chest through leadership, education, research, and communication.

This publication's content contains general information, is not intended to be and is not complete, is not medical advice, and does not replace professional medical care and physician advice, which always should be sought for any specific condition. The American College of Chest Physicians and its officers, regents, executive committee, members, and employees specifically disclaim all responsibility for any liability, damages (actual or consequential), loss, or risk, personal or otherwise, based on any legal theory whatsoever, alleged to have been incurred as a result, directly or indirectly, of the use of any of the material herein.