

# Philips Recall

Updated 6-28-21

## CPAPs, Bilevel PAPs, & Ventilators



### What Should I Do?

- Visit [philips.com/src-update](https://philips.com/src-update) to register your device for repairs or replacement.
- Stop using ozone-related cleaning products.
- Contact your physician.



### What Are the Health Risks?

- Philips has noted a breakdown of a type of foam in the device, resulting in breathing in foam particles and fumes.
- Side effects include headache; irritation of skin, eyes, and respiratory tract; effects on organs, such as kidneys and liver; as well as toxic and carcinogenic effects.



### Will My Insurance Cover a New Device?

- If your device is >5 years old, your insurance may pay for a new one. Please contact your physician.
- Otherwise, register the device with Philips for replacement or repairs.



### Should I Stop Using my Device?

- Philips recommends you stop using your CPAP/bilevel PAP.
- Contact your physician immediately to discuss your treatment options. Sudden discontinuation may have its own side effects.
- **Do not discontinue use of home ventilator devices**—contact your physician immediately for instructions. Use an inline bacterial filter.

For additional details, a full list of affected products, and potential health risks, visit [philips.com/src-update](https://philips.com/src-update) or call **877-907-7508**.